<b>BAUSCH &amp; LOMB</b> <b>IOM spa</b> <b>macherio (MI)</b>		PRODUCT					COUNTRY		<b>407439-051 DIS</b>	<b>APPROVAL SIGNATURE</b>	
		<b>RENU MULTIPLUS</b>					<b>RUS/UKR</b>				
<b>ARTWORK</b> <b>DEPARTMENT</b>		ARTWORK	REVISION	DCO N.	DRAWING SPEC.	DIMENSIONS	COMPONENT DESCRIPTION				
		<b>184082171</b>	<b>05</b>	<b>876S</b>	<b>55147</b>	<b>120 x 170</b>	<b>INSERT</b>				
		EUROLASER REFERENCE	USE WITH DESIGN ARTWORK N.			NOTE:	COLORI DI STAMPA TOTALI				
							<b>1</b>				
COLOUR 1	COLOUR 2	COLOUR 3	COLOUR 4	COLOUR 5	COLOUR 6	COLOUR 7					
<b>NERO</b>											

**GDP srl • TEL. 02.38.106.1**

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## ReNu MultiPlus® Multi-Purpose Solution

This insert contains important information. Please read it carefully and retain for future reference.

Bausch & Lomb ReNu MultiPlus Multi-Purpose Solution

**With HYDRANATE® Protein Remover**

Fresh Lens Comfort

Moisturizes, Cleans, Disinfects, Lubricates, Rinses, Stores and Removes Protein.

**CONTENTS:**

A sterile, isotonic solution that contains boric acid, edetate disodium, sodium borate and sodium chloride; DYMED® (polyamino-propyl biguanide) 0.0001%, HYDRANATE® (hydroxyalkylphosphonate) 0.03%, poloxamer 1%.

**ACTIONS:** Cleans, loosens and removes accumulations of film, protein and other deposits from soft contact lenses. DESTroys harmful micro-organisms on the surface of the lens. Rinses, stores and lubricates.

**INDICATIONS:** BAUSCH & LOMB ReNu MultiPlus® Multi-Purpose Solution is indicated for use in chemical (not heat) disinfection, daily cleaning, daily protein removal, lubricating, rinsing and storage of all soft contact lenses as recommended by your eye care practitioner.

**DIRECTIONS:**

Always wash and rinse your hands before you handle your lenses. Lenses must be thoroughly cleaned, rinsed and soaked each time they are removed to achieve disinfection.

**To clean, disinfect and remove protein daily, complete these simple steps. This daily regimen is recommended by Bausch & Lomb for a healthy and comfortable contact lens wearing experience:**

**Step 1:** Place 3 drops of ReNu MultiPlus Multi-Purpose Solution on each side of the lens surface and gently rub for 10 seconds on each side.

**Step 2:** Thoroughly rinse each side of the lens for 5 seconds with ReNu MultiPlus Multi-Purpose Solution.

**Step 3:** Place cleaned contact lens in the lens case and fill with fresh ReNu MultiPlus Multi-Purpose Solution. Soak at least 4 hours. Remember to always use fresh solution – never reuse solution.

Always follow your eye care professional's instructions. Based upon your individual tear chemistry and lens-wearing schedule, your eye care professional may recommend additional products or procedures. Your eye care professional may also recommend an alternative procedure without a rubbing step. If so, to clean, disinfect and remove protein daily without a rub step, complete Steps 2 and Step 3 only, as described above.

If any debris remains on your lenses you may rinse with ReNu MultiPlus Multi-Purpose Solution. Pour the BAUSCH & LOMB ReNu MultiPlus Multi-Purpose Solution out of the lens case and rinse the case with fresh ReNu MultiPlus Multi-Purpose Solution. Leave the case with the caps off to dry in the air. Replace the caps. Your case is now ready for use.

**STORAGE:** You may store your lenses in the unopened lens case until ready to wear, up to a maximum of 30 days. If you store your lenses for longer periods of time, they must be cleaned and disinfected with fresh ReNu MultiPlus Multi-Purpose Solution every 30 days and prior to wearing.

**Weekly Care**

No separate protein remover is required to remove protein deposits from the lenses. Daily use of BAUSCH & LOMB ReNu MultiPlus Multi-Purpose Solution as directed removes protein.

**WARNINGS & PRECAUTIONS:**

Lens care procedures recommended by your practitioner must be followed. Failure to follow these procedures may result in the development of serious eye infections.

- Always wash your hands before handling your lenses.
- Discard the BAUSCH & LOMB ReNu MultiPlus Multi-Purpose Solution from the lens case after each use.
- If you are allergic to any ingredient in this product, do not use.
- Discard solution six months after opening.
- Use before expiration date marked on the carton and bottle.
- Do not touch the bottle tip to any surface or your eye since this may contaminate the solution or cause injury to the eye. Replace the cap after use.
- Never re-use this solution.
- Do not use any eye medication in conjunction with this solution unless under medical supervision.
- Do not use BAUSCH & LOMB ReNu MultiPlus Multi-Purpose Solution with a heat disinfection method.
- Keep out of the reach of children.

**Important:** If irritation or excessive tearing occurs, persists or increases, or if vision is impaired, discontinue use and promptly consult your eye care practitioner.

184082171-05

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© Bausch & Lomb Incorporated.
BAUSCH & LOMB, ReNu MultiPlus, DYMED и/та and HYDRANATE являются зарегистрированными торговыми марками корпорации Бауш энд Ломб (Bausch & Lomb Inc.)/торгові марки Bausch & Lomb Incorporated/are trademarks of Bausch & Lomb Incorporated.